Stifneck Select

Select and Pedi-Select Extrication Collar

Directions for Use

980010 Stifneck Select (Qty:1)
980012 Stifneck Select, green (Qty:1)
980020 Stifneck Pedi-Select (Qty:1)
989910 Stifneck Select, Private Label (Qty: 1)
989920 Stifneck Pedi-Select, Private Label (Qty: 1)

Manufactured in China for:
LAERDAL MEDICAL AS, Norway
P.O. Box 377, N-4002 Stavanger, Norway
Tel. +47 51 51 17 00,
Fax +47 51 52 35 57

Distributed by:
Canada and Latin America:
LAERDAL MEDICAL CANADA LTD.
151 Nashdene Rd., Unit #45
Toronto, ON, Canada, M1V 4C3
Tel. +1 (416) 298-9600,
Toll free 888/LAERDAL (523-7325)
ou en français (800) 567-9987
Fax +1 (416) 298-8016

USA:
LAERDAL MEDICAL CORPORATION
167 Myers Corners Road, P.O. Box 1840
Wappingers Falls, New York 12590-8840
Tel. (800) 431-1055, +1 (845) 297-7770
Fax (800) 227-1143, +1 (845) 298-4545

Australia
Laerdal Pty. Ltd.
8 Stamford Road
Oakleigh VIC 3166.

Laerdal do Brasil
Al. Tocantins 125 – sl 01 – Alphaville
Barueri / SP / Brazil
06543-260
Tel: +55 11 4193-8007
e-mail: comercial@laerdal.com

www.laerdal.com

ENGLISH
ČESKÝ
DANSK
DEUTSCH
ESPANÖL
FRANÇAIS
ITALIANO
NEDERLANDS
NORSK
POLSKI
PORTUGUÊS
Русский
SVENSK
SUOMI
日本語
中文
한국어

Laerdal

helping save lives


Contents

| | |
|------------------------|----|
| ENGLISH | 4 |
| ČESKÝ | 5 |
| DANSK | 6 |
| DEUTSCH | 7 |
| ESPAÑOL | 8 |
| FRANÇAIS | 9 |
| ITALIANO | 10 |
| NEDERLANDS | 11 |
| NORSK | 12 |
| POLSKI | 13 |
| PORTUGUÊS | 14 |
| Русский | 15 |
| SVENSK | 16 |
| SUOMI | 17 |
| 日本語 | 18 |
| 中文 | 19 |
| 한국어 | 20 |
| REGULATORY INFORMATION | 21 |

Intended Use

This product is intended to support cervical spine in a neutral position during transportation, in combination with other cervical and full body immobilization devices (CSIDs and FBSIDs).

Cannula hooks are provided for the securement of supplemental oxygen lines.

Important


- In cases of a suspected spinal cord injury, proper neck immobilization is only one part of total immobilization.
- It is imperative that the patient be properly immobilized to prevent any movement of the spinal column.

Warnings

- Stifneck Select Collars should only be used by qualified healthcare providers.
- Improper use may result in injury or permanent disability. Always follow local spinal immobilization protocol.
- Use the correct size collar: An oversized collar may hyperextend a patient's cervical spine. An undersized collar may not provide appropriate stability.
- Do not adjust the Select or Pedi-Select collar on patient.

User Instructions


1. Measure the Patient


Align the head to neutral or "eyes forward" position unless contraindicated by local protocol.

2. Match the collar size to the patient:


Select: Choose from 4 adult positions.
Pedi-Select: Choose from 3 Paediatric positions.


3. Adjust and Lock the adjustable collar

Use the same method for adjusting both Collars.


4. Preform the collar


5. Apply the collar while manually maintaining neutral head position

Place the chin support under the chin. If a different size is needed, remove, resize, and re-apply the collar.

Pull the back of the collar snug while holding the front in place, and then fasten.

For a supine patient, slide the rear panel behind the neck before placing the chin support.

Technical Specifications

Storage:
 Do not store collar in folded position.
 Store in flat position.

Storage Temperature Range:
 -34 °C to 52 °C
 (-30 °F to 125 °F)

Operating Temperature Range:
 -25 °C to 43 °C
 (-13 °F to 110 °F)

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